

**OPERATOR**  
Instructions

# **PERIMATIC GP OPERATOR MANUAL**



**BS EN ISO 9002**  
**N° FM 15114**  
**N° RS 27480**

**CE**

**SOFTWARE VERSION 4.00**  
**DOC REV 1.3**

**JENCONS**  
(SCIENTIFIC) LTD

**MANUFACTURER & DISTRIBUTOR  
OF LABORATORY SUPPLIES  
EQUIPMENT • DISPOSABLES • CONSUMABLES**

<i>Chapter</i>	<i>Section</i>	<i>Page</i>
<b>Introduction</b>	.....	i
	Installing the Tubing .....	i
	Interfacing with the <b>PERIMATIC GP</b> .....	ii
	Approvals .....	iii
	Quickstart for a pre-programmed <b>PERIMATIC GP</b> .....	iii
<b>Keypad Layout</b> .....	1	
<b>Initial Setup</b>	.....	2
	Selecting a language .....	2
	Clearing the memory.....	2
	Exit Setup Mode .....	2
<b>Setting the Operating Parameters</b>	.....	3
	Preparing for dispensing .....	3
	Setting the dispensing volume .....	3
	Setting the repeat & delay between dispensings.....	4
	Automatic & Manual Modes .....	5
	Setting the pumping rate .....	6
	Setting the tubing diameter .....	7
	Setting the rotor direction .....	8
	Rounding up or down for greater accuracy.....	8
	Priming & calibrating the system .....	9
<b>Dispensing</b>	.....	11
	Manual & automatic dispensing .....	11
	Dispensing in conjunction with the <b>PERIMATIC ROBOTIC MODULE</b> .....	12
	Manual dispensing using the foot switch.....	12
<b>Memory Facilities</b>	.....	13
	Storing a program in memory .....	13
	Recalling a program from memory .....	14
	Reviewing programs in memory.....	14
<b>Error Messages</b> .....	15	
<b>Specification</b> .....	16	

## **INTRODUCTION**

**OPERATOR**

Instructions


The **JENCONS PERIMATIC GP** has been carefully designed and manufactured as a programmable automatic pump which can be used as a stand alone unit or interfaced with the **JENCONS PERIMATIC ROBOTIC MODULE** filling station.

**Packaging** This system is supplied complete with an instruction manual, mains lead, foot switch and a range of silicone tubing (3mm, 5mm & 8mm).

**Technical Data** A full technical specification is detailed on page 16.

**Installing the Tubing** **PERIMATIC** silicone tubing is fitted by loosening the red knob and lifting off the tubing clamp, both of which are illustrated in *figure 1*. The tubing is then laid into the arc of the rotor casting, facilitated by turning the rotor by hand as required, before clamping the tubing into place and tightening the red knob.


**NOTE.** Ensure that the clamp does not pinch or constrict the tubing through over-tightening of the red knob. The unit will not operate unless the door panel is down.



**Figure 1** Installing **PERIMATIC** Tubing

**Interfacing with the** The **PERIMATIC GP** has two sockets on the rear panel:-

- **PERIMATIC GP** Mains socket
- Foot switch/ **ROBOTIC MODULE** socket



**Figure 2** Rear Panel of **PERIMATIC GP**

- 1 Ensure that the voltage selector, found on the underside of the unit, is switched to the appropriate voltage, either 120V for 110-120V operation or 240V for 220-240V operation.
- 2 Plug the mains lead into mains socket at rear of unit.
- 3 Use the second socket for connecting the foot switch or the interface lead (for use with the **PERIMATIC ROBOTIC MODULE**).

**Additional Components Available**

Cat. No.	Component
256-065	<b>PERIMATIC ROBOTIC MODULE</b>
686-203	Tubing (bore x $\phi$ ) 3x7mm
686-210	Tubing (bore x $\phi$ ) 5x9mm
686-216	Tubing (bore x $\phi$ ) 8x12mm

**Approvals** The **PERIMATIC** pump conforms to the following requirements :-

- Machine safety regulations No 3073 (1992)
- Elec equipment for measurement/lab use BS EN 61010-1
- EMC directive BS EN50081-1 / BS EN50082-1 (1992)
- Low voltage directive (73/23/CEE)


**Quick-start for a pre-programmed PERIMATIC GP** In order to initiate a pump sequence on a pre-programmed **PERIMATIC GP**, the user should switch the unit on (switch **PERIMATIC GP** located on rear panel), and will be presented with the following screen display:

Software Version X.X
MEMORY SETUP PUMP

Having checked that the silicone tubing is correctly mounted around the rotor and is terminated at appropriate source and recipient vessels, press the 'PUMP' button or the <START> key to begin a fill cycle. The **PERIMATIC GP** will operate using the parameters last used or most recently set up in memory.

The remainder of the manual explains, in detail, how to program the system for a variety of filling applications using different tubing sizes and fluid types.

**Introduction** The front panel keypad is shown below in *figure 3*. Beneath the display are three rectangular soft function keys whose function is denoted by the text located directly above them, on the lower line of the LCD display.



**Figure 3** *Front Panel Keypad*

The keypad serves two functions:

- To enable parameters to be entered in the operating program of the **PERIMATIC GP**.
- To enable the operator to prime, calibrate and run the system.

**Conventions used in the manual** Throughout the manual the following conventions are used to describe the keyboard keys:

- The soft keys are referred to as ‘buttons’, e.g. from the main menu press the ‘MEMORY’ button.
- The keypad keys and the START/STOP keys are shown enclosed thus <START>.

**Initial Setup** This section illustrates how the user would set up various parameters, such as language used, speed etc, upon receipt of a new machine. It also shows that this procedure need only be carried out once as all parameters remain stored, even when the unit is switched off.

**Accessing the Setup Menu** The Setup Menu is accessed by holding down any one of the soft keys whilst the power is switched on at the rear of the machine. When the key is released the display will show the following message:

Setup Program  
Enter Password

Enter the system password (default 1212) the display will show the following message (Setup Menu):

Setup Program  
COUNTRY MEM ACCEPT

**Selecting a Language** Press the key beneath the 'COUNTRY' option on the display. The display will show the following message:

ENGLISH NEXT ACCEPT

Still using the soft keys, press 'NEXT' until the desired language is displayed and press 'ACCEPT' to confirm this selection. The program will then return to the Setup Menu.

**Clearing Memory** **NOTE.** This should only be necessary after maintenance work or after a software upgrade. This operation will erase all settings, including calibration, stored in the module.  
**NB.** Additional parameters can be set up by an Engineer for special applications. Clearing the memory will also clear these settings back to defaults.

At the setup menu press the 'MEM' button. The display will show the following message:

Erase Memory  
ACCEPT QUIT

Press the 'ACCEPT' button to erase all memory locations, or the 'QUIT' button to cancel the operation and return to the Setup Menu.

**Exit Program Mode** To exit the Setup Menu press the 'ACCEPT' button in the main Setup Menu above. The **PERIMATIC GP** will reset and the display will show the Main Menu.

## SETTING THE OPERATING PARAMETERS

**Introduction** Before dispensing is performed the **PERIMATIC GP** needs to know certain things, these are :-

- 1 How much liquid to dispense.
- 2 Whether the dispenser is running in Automatic or Manual Mode.
- 3 How many times to repeat dispensing (Automatic Mode only).
- 4 The delay from the finish of one dispense to the start of the next dispense (Automatic Mode only).
- 5 The pumping profile and timings.
- 6 The tubing diameter in use.
- 7 Whether better accuracy is required by pumping an exact number of rotor half turns.
- 8 The direction of rotation.

The **PERIMATIC GP** is switched on by depressing the mains ON/OFF switch, located above the mains socket on the rear panel.

When the **PERIMATIC GP** has initialised the display will show:

Software Version X.X
MEMORY    SETUP    PUMP

X.X identifies the installed software version number.

Press 'SETUP'; the display will show:

Setup GP
PRIME            CAL    OPTIONS

Press 'OPTIONS'; the display will show:

Setup GP
VOLUME    COUNT    NEXT

**Setting the dispensing volume** Press 'VOLUME'; the display will show:

Volume	10.0 ml
ACCEPT    QUIT    DELETE	

The previously set volume will be displayed, enter the volume of fluid required to be dispensed (in millilitres) using 'DELETE' to erase any erroneous entries.

Press <ENTER> to update the display with the new value.

Press 'ACCEPT' to store the volume and exit to the previous menu.

Press 'QUIT' to exit from this menu without changing the original volume.

**NOTE.** Using 'QUIT' is a quick way of verifying the programmed volume.

**Setting the repeat dispense count** From the Setup Options Menu press 'COUNT'; the display will show:

No of Doses 10
ACCEPT    QUIT    DELETE

The previously set value will be displayed. Enter the number of times that you require the dispensing to be repeated when Automatic Mode is in use. Use 'DELETE' to correct any errors.

Press <ENTER> to update the display with your entered value.

Press 'ACCEPT' to store the value and exit to the previous menu or 'QUIT' to exit from this menu without changing the original value.

From the Setup Options Menu press 'NEXT'; the display will show:

Setup GP
DELAY    MODE    NEXT

**Setting the delay between dispensings** When the **PERIMATIC GP** is operating in the Automatic Mode a delay is inserted between dispensings. This delay can be set from 0.1 to 9.9 seconds.

Press 'DELAY'; the display will show:

Dose Delay 1.0 s
ACCEPT    QUIT

Enter the time delay using the cursor keys (the decimal point is inserted automatically).

Press 'ACCEPT' to exit and store the value or 'QUIT' to exit without changing anything. The previous menu will be displayed.

**Automatic & Manual Modes** In Automatic Mode the dispenser will perform a dispense (set using 'VOLUME' command) and then wait for a predefined time (set using the 'DELAY' command). This sequence may then be repeated a number of times (set using 'COUNT').

In Manual Mode the **PERIMATIC GP** will dispense one dose and return to its previous menu.

**Selecting Automatic or Manual Modes** Press 'MODE'; the display will show:

```
Dispensing is manual
NEXT    ACCEPT
```

The current selected mode (in this case Manual) is displayed on the top line.

Press 'NEXT' to change to the other (Automatic) mode; the display will show:

```
Dispensing is auto
NEXT    ACCEPT
```

Press 'ACCEPT' to exit and store the selected mode. The previous menu will be displayed.

Press 'NEXT'; the display will show:

```
Setup GP
PROFILE  TUBE  NEXT
```

**Setting the pumping rate** The pumping rate initially takes two forms; slow pumping and fast pumping. In addition to this the **PERIMATIC GP** also performs the following combinations :

START	RUN	FINISH
SLOW	SLOW	SLOW
FAST	FAST	FAST
SLOW	FAST	FAST
SLOW	FAST	SLOW
FAST	FAST	SLOW

**NOTE.** Using a slow start & finish reduces the possibility of splashing.

Press 'PROFILE'; the display will show:

<b>SLOW</b>	<b>SLOW</b>	<b>SLOW</b>
	<b>NEXT</b>	<b>ACCEPT</b>

The word on the left of the top line denotes the start speed, the middle word the run speed and the right word the finish speed.

Press the 'NEXT' button until the required combination of start, run & finish speeds are displayed.

You will notice that when a combination of slow and fast speeds are selected an additional 'TIME' button appears. This allows the time spent at the selected start or finish speed to be set.

**EXAMPLE.** Using the SLOW-FAST-SLOW speed option.

Press 'TIME'; the display will show:

<b>Start/ Finish Time</b>		
<b>START</b>	<b>QUIT</b>	<b>FINISH</b>

'START' allows the time spent at the selected start speed to be set. 'FINISH' allows the time spent at the selected finish speed to be set.

Press 'START'; the display will show:

<b>Start Time 1.0 s</b> <b>ACCEPT QUIT</b>
---

Enter the time to spend at the selected start speed between 0.1 and 9.8 seconds.

Press 'ACCEPT' to store the value and exit.

Press 'QUIT' to exit without changing the previous value.

The previous menu will be displayed.

Press 'FINISH'; the display will show:

<b>Finish Time 1.0 s</b> <b>ACCEPT QUIT</b>
--

Enter the time to spend at the selected finish speed; between 0.1 - 9.8 seconds.

Press 'ACCEPT' to store the value and exit to the previous menu or QUIT to exit without changing the previous value.

Press 'QUIT'; the display will show:

<b>SLOW</b>	<b>FAST</b>	<b>SLOW</b>
<b>TIME</b>	<b>NEXT</b>	<b>ACCEPT</b>

Press 'ACCEPT' to store the settings. The previous menu will be displayed.

**Setting the tubing diameter** To enable the **PERIMATIC GP** to dispense fluid accurately it needs to know the size of the tubing in use. These sizes are specified by the **internal diameter** (i.d.) of the tube. The **PERIMATIC GP** can accept tubing in sizes from i.d. 3.0 to 8.0mm.

Press 'TUBE'; the display will show:

<b>3mm I.D.</b>	<b>NEXT</b>	<b>ACCEPT</b>
-----------------	-------------	---------------

The currently selected tube diameter is displayed. Press 'NEXT' until the correct tubing size is displayed.

Press 'ACCEPT' to store the selection and exit to the previous menu.

Press 'NEXT'; the display will show:

Setup GP  
DIRECTION ROUND NEXT

**Setting the rotor direction** The rotor on the **PERIMATIC GP** can run either clockwise or anti-clockwise.

Press 'DIRECTION'; the display will show:

Run Clockwise  
NEXT ACCEPT

Press 'NEXT' to change the rotation direction and 'ACCEPT' to store the setting. The previous menu will be displayed.

**Rounding up or down for greater accuracy** The **PERIMATIC GP** will dispense most volumes with good repeatability. However, due to the fact that the start and finish points may not always fall on a complete half turn of the rotor, errors may be apparent.

The 'ROUND' command allows the user to set the rotor movement to a multiple of complete half turns. For example. If a pump is dispensing and the rotor takes  $2\frac{1}{4}$  turns rounding up will make the rotor take  $2\frac{1}{2}$  turns and rounding down will take 2 turns.

Press 'ROUND'; the display will show:

Round to Half Turn  
↓↓ QUIT ↑↑

Press UP '↑↑' to round the volume up to the nearest half rotor turn and DOWN '↓↓' to round down to the nearest half rotor turn.

**NOTE.** When the **PERIMATIC GP** is dispensing, the display will show 'Rounded' if the round option has been selected. The display will reflect the calculated, rounded volume as the dispense volume.

This completes the set up of **PERIMATIC GP** options. Press <ENTER> to return to the setup menu. The display will show:

Setup GP		
PRIME	CAL	OPTIONS

**Priming the system** Press 'PRIME'; the display will show:

Ready to Prime
ACCEPT    QUIT

Once the 'ACCEPT' button is pressed the rotor will begin to turn and fluid will be drawn into the tubing. Press 'ACCEPT' when the tubing is full and the rotor will stop.

Press 'QUIT' to return to the main menu.

**Calibration procedures** The **PERIMATIC GP** has default settings for flow constants with various tubing sizes. These are based on measurements taken using water at 20°C. The **PERIMATIC GP** can be recalibrated to take account of liquids that are at a higher or lower temperature and/or have different viscosities.

**NOTE:** The calibration process requires an accurate balance to be available to ascertain the weight of the dispensed volume of liquid. Furthermore, calculations may be necessary to convert from weight in grams to millilitres based on the viscosity and temperature of the liquid.

**NOTE:** Modification to the calibration constants are made on the currently selected tubing diameter only. If different tubing sizes are used then each must be calibrated if required.

From the main menu press 'SETUP'; the display will show:

Setup GP		
PRIME	CAL	OPTIONS

Press 'CAL'; the display will show:

<b>Calibrate GP</b>
<b>VOLUME    QUIT    ROUND</b>

Press 'VOLUME'; the display will show:

<b>Volume    10.0 ml</b>
<b>ACCEPT    QUIT    DELETE</b>

The previously set calibration volume will be displayed. Enter the volume required for calibration, in millilitres. Use 'DELETE' to correct any errors.

Press <ENTER> to update the display or 'QUIT' to exit from this menu without changing the currently set volume.

Press 'ACCEPT' to store the volume and exit to the next menu as shown:

<b>Press START to pump</b>
<b>QUIT</b>

Press <START> to dispense the volume entered. When the **PERIMATIC GP** has finished pumping, the display will show:

<b>Actual Vol    10.0 ml</b>
<b>ACCEPT    QUIT    DELETE</b>

Ascertain the volume of liquid by weighing (allow for temperature & viscosity). Enter the actual volume dispensed in millilitres. Use 'DELETE' to correct any errors.

Press <ENTER> to update the display and 'ACCEPT' to store the volume and exit to the next menu.

The **PERIMATIC GP** will recalibrate its new flow constants for the liquid and then return the user to the 'PRESS START to PUMP' display.

Press 'QUIT' to exit from this menu and store the calibration values and, finally, press <START> to pump a new sample to check the dispensed volume.

**Manual dispensing** In Manual Mode press 'PUMP' or <START> to display:

1	10.0ml
MAN	SFS

The **PERIMATIC GP** will dispense its programmed volume, using the programmed speed profile.

During manual dispensing the display shows the following information :

- 1 Total number of manual dispensings since setup. Shown on the top, left side of the display.
- 2 The dispense volume. Shown on the right side of the top display line.
- 3 The dispense mode (MAN). Shown on the left side of the bottom display line.
- 4 The pumping rate; shown in the middle of the bottom display line.
- 5 'Rounded' will be displayed in the bottom right corner of the display if the 'ROUND' option has been selected.

When the **PERIMATIC GP** has finished dispensing the display reverts to the main menu. Subsequent pressing of either 'PUMP' or <START> will repeat the dispensing process.

**NOTE.** The dispense count will increment each time 'PUMP' or <START> is pressed.

**Automatic dispensing** In Automatic Mode press the 'PUMP' button or <START> key and the display will show:

1 of 100	10.0ml
AUTO	SFS

The **PERIMATIC GP** will dispense the programmed volume, using the programmed speed profile, and wait for the programmed delay time.

During automatic dispensing the display will show the following information :

- 1 A continuous tally count of dispenses completed and the programmed number of dispenses to reach is shown on the top left side of the display (e.g. 1 of 100).
- 2 The dispense volume is shown on the top, right side of the display.
- 3 The dispensing mode (AUTO); shown on the bottom, left side display.
- 4 The speed profile is shown in the middle of the bottom display line.
- 5 'Rounded' will be displayed in the bottom right corner of the display if the 'ROUND' option has been selected.

When the programmed number of dispenses are complete the display reverts to the main menu. Subsequent use of the 'PUMP' button or <START> key repeats the dispensing process.

- Dispensing with the PERIMATIC ROBOTIC MODULE**
- 1 Set up the **PERIMATIC GP** for manual dispensing.
  - 2 Set up the **PERIMATIC ROBOTIC MODULE** for the appropriate fill profile (see separate instruction manual).
  - 3 Connect the **PERIMATIC GP** to the **PERIMATIC ROBOTIC MODULE** using the interface lead.
  - 4 With the Main Menu shown on both units, press <START> on the **PERIMATIC ROBOTIC MODULE**.

The **PERIMATIC ROBOTIC MODULE** will move to its first start position and the **PERIMATIC GP** will dispense. The **PERIMATIC ROBOTIC MODULE** will then continue to control the **PERIMATIC GP** for all other programmed fill positions.

When the fill sequence is finished both units will stop and revert to their Main Menus.

- Dispensing using the foot switch**
- Plug the foot switch into the 3 pin socket on the rear of the **PERIMATIC GP**. Press the foot switch once to start the pumping operation and again to stop. Using the foot switch in this way is identical to using the 'PUMP' soft key or <START> button in Manual or Automatic Mode.

When power is removed from the **PERIMATIC GP** the current setup is remembered. In addition to this the **PERIMATIC GP** can store and recall up to 10 previous settings, each of which may have a 20 character description attached to it to describe its various parameters.

All operating information for the **PERIMATIC GP** is stored in a memory location, including the language in use as well as pump speed and tubing flow constants.

From the Main Menu press 'MEMORY' to enter the Memory Menu:

Recall/ Store/ Review  
RECALL STORE REVIEW

**Storing a program in memory** From the Memory Menu press 'STORE':

Store Memory No 01  
ACCEPT QUIT

Enter the memory number to which you wish to store the current setup using the numeric keypad.

Press 'QUIT' to return to the main menu. Press 'ACCEPT' to store the program then the memory description screen will appear:

01:PROGRAM TITLE  
ACCEPT QUIT

The display will show the current description stored for that location. Use keys <4> and <6> to move cursor left and right and the <2> and <8> keys to increase or decrease the character under the cursor. Pressing <5> will copy the character to the left of the cursor to the current cursor position.

Press 'ACCEPT' to store the description and return to the main menu.

**Recalling a program from memory** From the Memory Menu press 'RECALL':

Recall Memory No 01  
ACCEPT QUIT

Enter the memory number to recall using the numeric keypad. Press 'ACCEPT' to recall the program or 'QUIT' to return to the Main Menu.

If the memory number entered does not contain a valid program the display will show:

\*\*\*\*\* ERROR \*\*\*\*\*  
MEMORY IS EMPTY

Re-enter a valid memory location number or use 'REVIEW' (below) to select from the valid memories stored.

**Reviewing programs in memory** From the Memory Menu press 'REVIEW':

01:PROGRAM TITLE  
NEXT QUIT RECALL

This will show the description attached to the first valid memory location. Pressing 'NEXT' moves to the next location with a valid program stored in it. When the last stored location is displayed and 'NEXT' is pressed the first location is displayed again.

If 'REVIEW' is pressed and there are no valid programs in memory then the review button appears to be ignored.

Press 'RECALL' to load the contents of the selected memory location (you will be returned to the Main Menu) or press 'QUIT' to return to the previous menu.

Error Message	Description
*****ERROR***** MEMORY IS EMPTY	An empty memory location was selected. Re-enter a valid memory number or use re-view to select.
*****ERROR***** ONLY 01 - 10 ALLOWED	A memory number outside of the range 01 to 10 was entered (invalid). Re-type with a number between 01 and 10.
*****ERROR***** Cannot Pump 0ml	A dispensing volume of 0ml was entered (invalid). Try again with a volume between 0.5ml and 9999.9ml.
*****ERROR***** MAX VOL is 9999.9ml	The <b>PERIMATIC GP</b> can only dispense up to 9999.9ml. Re-enter a lower figure.
*****ERROR***** ONLY 0.1 - 9.9 ALLOWED	A time delay was entered as 0 (invalid). Re-type with a number between 0.1 and 9.9.
*****ERROR***** Max doses are 999	The <b>PERIMATIC GP</b> can only repeat a dispense count of up to 999. Re-enter a lower figure.
Operating Error	Description
Speed profile set to SLOW /FAST/ SLOW but <b>PERIMATIC GP</b> pumps all of dose on fast	SPEED profile has been selected but not the START and FINISH times.
Speed profile set to SLOW /FAST/ SLOW but <b>PERIMATIC GP</b> pumps runs at slow all of the time.	Time allocated for the START SLOW run is too long OR the volume specified is dispensed before the slow run time has expired.

## ***SPECIFICATION***

***OPERATOR***  
*Instructions*

**Accuracy** -  $\pm 1.0\%$  for volumes over 20ml

**Repeatability** - Generally  $<1\%$  Discrepancy

**Dispensing volume** - 0.5ml (min) 9999.9ml (max)

**Tubing bore** - 3mm (min) 8mm (max)

**Tubing wall** - 2mm thick (max)

**Pumping method** - Peristaltic

**Rotor speed** - 70 rpm (slow)  
- 170 rpm (fast) for pipe size 3mm to 5mm  
- 150 rpm (fast) for pipe size 6mm to 8mm

### **Pumping rate**

**3mm bore** - 2.1ml/s (slow) 4.4ml/s (fast)  
**5mm bore** - 4.6ml/s (slow) 11.8ml/s (fast)  
**8mm bore** - 11.9ml/s (slow) 22.4ml/s (fast)

**Electronics** - Microprocessor controlled

**Program Storage** - 10

**Languages** - English, French, German, Italian, & Spanish  
(selectable)

**Power supply** - 110/120 and 220/240V (selectable)  
50/60Hz single phase

**Power consumption** - 90W


**Fuse rating** - T1.6A

**I/O ports** - Foot switch/XY Module

**Dimensions** - 290 x 250 x 150mm (WxDxH)

**Weight** - 7Kg (net)

## **PERIMATIC GP**



**PERIMATIC ROBOTIC MODULE** - available as option

**JENCONS**  
(SCIENTIFIC) LTD

☐ **Jencons (Scientific) Limited**  
Cherrycourt Way Industrial  
Estate, Stanbridge Road,  
Leighton Buzzard, Bedfordshire,  
LU7 4UA, ENGLAND

✉ techsales@jencons.co.uk

☎ +44 (0)1525 372010

📄 +44 (0)1525 379 547

☐ **Jencons Scientific Inc.**  
Bursca Business Park,  
800 Bursca Drive,  
Suite 801,  
Bridgeville PA 15017, USA

✉ info@jencons.com

☎ 1-412-257-8861

800 846-9959 (Toll Free)

📄 1-412-257-8809